

LEAD MAGNET | BOULDER BIOLABS

EO Residuals Sample Submission Guide

For ISO 10993-7 compliance, sample selection, and submission.

Co-located with Boulder Sterilization. One ISO 13485 quality system. No shipping samples across the country.

EO Residuals Sample Submission Guide

How to sample, label, and ship for ISO 10993-7 compliance - or just hand them across the campus.

Sample Selection

Patient contact category determines the limit.

Limited (<24h), prolonged (24h-30d), permanent (>30d) - each has its own threshold. Confirm your category before submitting.

Sample size guidance.

Typical: 3-10 devices per extraction condition, depending on validation phase. Initial validation needs more, annual revalidation often needs fewer.

Extraction strategy.

Headspace simulates patient use; exhaustive recovers total residual load. Submit guidance per ISO 10993-7 Section B.

Hold time matters.

Sample at end of aeration. Every additional hour of out-gassing pre-test inflates your validated aeration time.

Submission Checklist

- Device sterilized (and aerated) on documented date/time
- Lot number recorded
- Cycle parameters attached (aeration time, temperature, EO concentration)
- Patient contact category specified
- Validated ISO 10993-7 limits noted
- Samples packed in EO-permeable wrap or sealed pouches per device protocol
- Chain of custody attached (sterilization lot release record)
- Submission form completed

For Boulder Sterilization clients: the chain-of-custody and aeration documentation transfers internally - no chain-of-custody paperwork required from your team.

How to Reach Us

Email submissions: info@boulderiq.com

Schedule a pickup: [303-531-1238](tel:303-531-1238)

Quote requests: boulderbiolabs.com/contact